

ART-DIS-01 Disinfection guide card

PKG-EEG-015 Rev B
generated by tools/artwork_gen.py
2026-09-02 CC BY-SA 4.0

Size A5 portrait, 148 x 210 mm, laminated
Placement CONSUMABLES bay, packing-list line 6.6 (RFQ A-05)
Source SVC-EEG-013 Rev B section 2 R5, which section 9 of that document names as this card's source text

CLEANING AND DISINFECTION

Operator card. Extract of SVC-EEG-013 Rev B section 2 R5.
Efficacy is not validated by this programme (SVC-EEG-013 section 2 R5).

Item	Agent	Concentration	Contact	Method
8 cups, 2 ear clips	Neutral enzymatic instrument detergent	0.5 % v/v in demineralised water	10 min	Ultrasonic, 40 kHz
8 cups, 2 ear clips	Isopropanol	70 % v/v	2 min wet	Immersion
Boom capsule	Isopropanol	70 % v/v	2 min wet	Wipe only, never immersed
Boom windscreen	--	--	--	Discard and replace
HM-01 frame, HM-03 occipital yoke, HM-06 chin strap	Isopropanol	70 % v/v	2 min wet	Wipe, never immersed
Case interior and exterior	Quaternary-ammonium as supplied surface wipe or 70 % IPA		per label, min 1 min	Wipe
Headphone ear cushions	70 % IPA	70 % v/v	1 min	Wipe

NEVER

- Autoclave, any part - MJF PA12 and TPU 85A distort; the +0.15 / -0.05 mm bayonet fit is lost
- Immersion of HM-01 - The frame contains the harness, the LEDs and, in Phase 2, the electronics
- Ultrasonic on the boom capsule - The electret capsule is not sealed
- Bleach or any hypochlorite on Ag/AgCl - Destroys the sintered surface, and with it the lead-off calibration T10 depends on
- Acetone, MEK, any ketone - Crazes PA12 and attacks the PTFE-jacketed harness terminations
- Abrasives or wire wool on a cup face - Removes sintered material; the cup is no longer sintered, it is plated
- Dishwasher - Alkaline detergent attacks Ag/AgCl and the TPU

Record the agent, its lot and its expiry on F1 section 4. An expired agent invalidates the turnaround and the cleaning is repeated.

ART-DIS-01 PKG-EEG-015 Rev B source SVC-EEG-013 Rev B CC BY-SA 4.0

trim 148 x 210 mm, shown 1:1

Notes

- The agent, concentration and contact-time table and the prohibition table are read out of SVC-EEG-013 at generation time. Nothing on this card is written here, so the card and the manual cannot disagree.